




Disposable Suction & Irrigation Unit



Disposable L Hook & Spatula
for suction irrigation unit



Specifications

SI06-360T $\phi 5 \times 360$ mm

SI06-430T $\phi 5 \times 430$ mm

Product Features

- ◆ Ergonomically designed handpiece can be easily operated by either hand
- ◆ Unblocked channel, good flowing performance
- ◆ Stainless steel probe with position function, easy to fix and draw out
- ◆ Flip designed rear port for easy access inside chamber, can be operated by one hand, quick close and reduce air leakage
- ◆ Transparent headpiece, unobstructed visualization
- ◆ Resistant structure and reliable quality
- ◆ Clear channel, consistent diameter with self-flush function
- ◆ Smoke evacuation and filtration, smoke flow rate is adjustable

Specifications

LH07-360 $\phi 5 \times 360$ mm

Sp07 -360 $\phi 5 \times 360$ mm