


Polymer Ligating Clips and Appliers



Clips

Size	Colour	Vessel Range	Pcs per cartridge	Cartridges per box
Medium-Large Polymer Clips	Green	3 mm - 10 mm	6	14
Large Polymer Clips	Purple	5 mm - 13 mm	6	14
Extra-Large Polymer Clips	Yellow	7 mm - 16 mm	6	14



Appliers (Reusable)

Size

Open Applier Medium-Large	25cm
Open Applier Large	25cm
Open Applier Extra-Large	25cm

Endoscopic Applier Medium-Large	
Endoscopic Applier Large	
Endoscopic Applier Extra-Large	

Remover Medium-Large	
Remover Large	
Remover Extra-Large	

+ Surgical Specialities

- * General
- * Colorectal
- * Gynaecology

Welfare Polymer Clips are made from a permanent, non-absorbable polymer that is inert, non-toxic and radiolucent.